


## mRNA Capture Kit

for capturing 192 poly(A)<sup>+</sup> RNA preparations in PCR tubes


Cat. No. 11 787 896 001

<b>Principle</b>	<p>A lysis step releases RNA from cultured cells or tissues. The poly(A)<sup>+</sup> tail of mRNA in the lysate hybridizes to biotin-labeled oligo(dT). A streptavidin-coated PCR tube immobilizes the biotin-labeled dT-A hybrids and washes remove unbound contaminants. The oligo(dT) part of this hybrid also serves as a primer for reverse transcriptase in RT-PCR. The entire process (including RT-PCR) takes place in a single tube.</p> <p> <i>The mRNA Capture Kit contains only reagents for the purification of mRNA; it includes neither reverse transcriptase nor reagents for amplification.</i></p>
<b>Starting material</b>	<p>Research samples may contain up to:</p> <ul style="list-style-type: none"> <li>● 40 µg total RNA</li> <li>● 5 x 10<sup>5</sup> cultured cells</li> <li>● 20 mg tissue</li> </ul>
<b>Application</b>	Preparation of highly purified poly(A) <sup>+</sup> RNA, which may be used directly for qualitative or quantitative RT-PCR
<b>Time required</b>	<ul style="list-style-type: none"> <li>● Total time: approx. 30 min (starting from total RNA)</li> <li>● Hands-on time: &lt;25 min (starting from total RNA)</li> </ul>
<b>Results</b>	<ul style="list-style-type: none"> <li>● Yield: Variable, depending upon starting material (See Part IV of “How to use the kit” in this article.)</li> <li>● Purity: When isolated mRNA is used as template for RT-PCR, the amplicons are clearly visible and free of background (See “Typical results with the kit” in this article.)</li> </ul>
<b>Benefits</b>	<ul style="list-style-type: none"> <li>● <b>Significantly reduces the time needed to isolate mRNA</b>, because the kit uses a semi-automated process that easily handles multiple samples</li> <li>● <b>Minimizes sample loss, handling time, and contamination</b>, because the kit allows rapid, efficient isolation of mRNA and RT-PCR in a single tube</li> <li>● <b>Improves the reliability and reproducibility of RT-PCR</b>, because the kit produces excellent RT-PCR templates, even from small amounts of starting material or material containing low-abundance mRNAs</li> </ul>


# 5

## How to use the kit

### I. Flow diagram



### II. Kit contents

- Lysis Buffer containing lithium dodecyl sulfate (50 ml)
- Oligo(dT)<sub>20</sub> Concentrate: Oligo(dT)<sub>20</sub>, biotin-labeled, 20x concentrated (50 µl)
-  Dilute the Oligo(dT)<sub>20</sub> Concentrate just before use.
- Double dist. Water, nuclease-free (1 ml)
- Wash Buffer, ready-to-use (500 ml)
- Streptavidin-coated PCR Tubes (192)
- Caps for PCR Tubes (192)

5

### III. Additional materials needed

- Sterile cups (for preparing dilutions)
- Pipettes with sterile, aerosol-resistant pipette tips
- Syringe fitted with a 21-gauge needle (for tissue and cultured cells only)
- PBS (for preparing cultured cells only)
- Mortar, pestle, liquid nitrogen (for preparing tissue samples only)
- PCR thermal cycler
- Reagents for RT-PCR
- Electrophoresis equipment or other equipment for analyzing PCR products

### IV. Typical mRNA content of different samples


Starting material	Total RNA (µg)	mRNA (µg)
Cultured cells (10 <sup>7</sup> cells)	30 – 500	0.3 – 25
Tissue (100 mg):		
Mouse brain	200	7
Mouse liver	700	14
Mouse lung	130	10

### V. Protocol for capturing mRNA from total RNA, cultured cells, or tissue

#### General notes:


- ▶ Make sure that all material which contacts the RNA is free of contaminating RNases. For decontamination of equipment, see Appendix (page 219).
- ▶ For tips on preventing RNase contamination in this procedure, as well as guidelines for the subsequent RT-PCR procedure, see the package insert provided with the mRNA Capture Kit.

- 1 To prepare the Oligo(dT)<sub>20</sub> Working Solution used in Step 3 below, dilute the Oligo(dT)<sub>20</sub> Concentrate 1:20 with nuclease-free double dist. water.

 Prepare the Oligo(dT)<sub>20</sub> Working Solution just before use; prepare only as much as you will need. (Each capture requires 4 µl of Working Solution.)

- 2 Prepare the sample:

- ▶ Total RNA: Dilute ≤40 µg of total RNA in 200 µl of Lysis Buffer.
- ▶ Cultured cells: Wash cells (≤5 x 10<sup>6</sup>) twice with ice-cold phosphate-buffered saline (PBS). Add 200 µl of Lysis Buffer to the cell pellet. Shear DNA mechanically by passing the sample 6 x through a 21-gauge needle.
- ▶ Tissue: Snap freeze ≤20 mg of tissue. Grind frozen tissue to a homogeneous powder in a pre-cooled mortar. Chill 200 µl of Lysis Buffer to ≤0°C in a sodium chloride-ice water bath, then add the frozen powder to the chilled Lysis Buffer. Homogenize the powder suspension by passing it 4 x through a 21-gauge needle. Centrifuge the suspension at 11,000 x g for 30 s. Use only the supernatant for the capture procedure.

 If possible, perform all steps in the capture procedure below at a temperature between 0°C and -4°C. Otherwise, perform these steps at 4°C.




- 
- 3** Hybridize the mRNA to biotin-labeled oligo(dT)<sub>20</sub>:
- ▶ Add 4 µl Oligo(dT)<sub>20</sub> Working Solution to the sample to form the hybridization mix.
  - ▶ Incubate the hybridization mix for:
    - ▶ EITHER 5 min at 2 to 8°C (if you started with tissue)
    - ▶ OR 5 min at 37°C (if you started with cultured cells or total RNA)
- 
- 4** Immobilize mRNA in the Streptavidin-coated PCR Tubes:
- ▶ Add 50 µl of the hybridization mix to a Streptavidin-coated PCR tube.
    - ! Use 4 tubes to process entire hybridization mix.
  - ▶ Incubate the hybridization mix for:
    - ▶ EITHER 3 min at 2 to 8°C (if you started with tissue)
    - ▶ OR 3 min at 37°C (if you started with cultured cells or total RNA)
- 
- 5** Wash the immobilized sample 3 times. For each wash:
- ▶ Add 250 µl Wash Buffer to the tube.
    - ! To avoid losing mRNA, do not wash the sample too vigorously.
  - ▶ Remove all Wash Buffer and discard.
    - ! After the third wash, the captured mRNA in the Streptavidin-coated PCR Tube is ready for RT-PCR.
- 

## VI. Troubleshooting the mRNA Capture protocol

If you get...	Then, the cause may be...	And you should...
No PCR product (after RT-PCR)	Insufficient homogenization of sample material	<ul style="list-style-type: none"> <li>▶ Eliminate all visible clumps in lysate (Step 1) by homogenization, etc.</li> <li>! mRNA trapped inside sample clumps can not be isolated by this procedure.</li> </ul>
	Too much starting material per volume of Lysis Buffer (lysate viscosity too high)	<ul style="list-style-type: none"> <li>▶ Do not use more than the recommended amount of starting material.</li> <li>▶ Dilute the lysate (Step 1) before continuing the procedure.</li> <li>! The viscosity of the lysate is critical to successful mRNA isolation.</li> </ul>
	Degradation of RNA by contaminating RNase from buffers or equipment	<ul style="list-style-type: none"> <li>▶ Check all buffers for RNase contamination.</li> <li>▶ Follow procedures in the package insert for RNase decontamination.</li> </ul>
	Degradation of RNA because RNases in sample material were not sufficiently inactivated	<ul style="list-style-type: none"> <li>▶ Follow the kit procedure exactly as written and be especially careful to:               <ul style="list-style-type: none"> <li>▶ Work rapidly</li> <li>▶ Perform indicated steps at ice water temperatures (0 to -4°C)</li> </ul> </li> </ul>

5

## Typical result with the kit



**Figure 48: Detection of MCAD and  $\beta$ -actin transcripts in K562 cells.** Lysates from  $10^6$  cells were serially diluted. Messenger RNA from the dilutions was captured in Streptavidin-coated PCR tubes according to the above protocol. RT-PCR was performed in the same PCR tubes with reagents from the Titan One Tube PCR System and primers derived from the human MCAD (Medium-chain acyl-CoA dehydrogenase) and  $\beta$ -actin genes. Amplicons were run on a 1% agarose gel. For details of the procedure and the primers used, see Zoelch and Frey (1996).

**Lanes 1, 10, 18:** DNA Molecular Weight Marker VI  
The remaining lanes contained mRNA from

**Lane 2:**  $1 \times 10^4$  cells

**Lanes 3, 11:**  $2 \times 10^3$  cells

**Lanes 4, 12:** 400 cells

**Lanes 5, 13:** 80 cells

**Lanes 6, 14:** 16 cells

**Lanes 7, 15:** 3.2 cells

**Lanes 8, 16:** 0.64 cells

**Lanes 9, 17:** Control (no RT performed)

**Result:** The mRNA Capture Kit and Titan One Tube RT-PCR System allow the efficient detection of  $\beta$ -actin and MCAD template RNA's from as few as 3.2 K562 cells.

## References

- Blümel, J. et al. (2005) *J. Virol.*, **79**, 14197 – 14206
- Breiner, H. W. et al. (2005) *Biochemica* **1**: 7 – 8
- Cloosen, S. et al. (2007) *Cancer Res.*, **67**, 3919 – 3926
- Cote, G. et al. (1997) *Biochemica* **4**, 14 – 18
- Feron, R. et al. (2004) *Biochemica* **3**: 23 – 24
- Fillingham, J. S. et al. (2006) *Eukaryot. Cell*, **5**, 1347 – 1359
- Fujita, M. et al. (2005) *Carcinogenesis*, **26**, 271 – 279
- García-Vallejo, J. J. et al. (2005) *J. Biol. Chem.*, **280**, 12676 – 12682
- Gringhuis, S. I. et al. (2005) *Mol. Cell. Biol.*, **25**, 6454 – 6463
- Jin, S. et al. (2005) *J. Virol.*, **79**, 8793 – 8801
- Vanderheyden, M. et al. (2005) *Heart*, **91**, 926 – 931
- Weitkamp, J.-H. et al. (2005) *J. Immunol.*, **174**, 3454 – 3460
- Zoelch, C. and Frey, B. (1996) *Application Note, The Journal of NIH Research*, No. 6