

### III. Buffers and Gels for Electrophoresis

#### Running buffers for non-denaturing gel electrophoresis

(Adapted from Brown, 1991)

Buffer	Components (1 x buffer)	Recipe for 1 liter of 10 x buffer
TBE <sup>a,b</sup> (Tris-borate)	89 mM Tris-borate (pH 8.3) 2 mM Na <sub>2</sub> EDTA	108.0 g Tris-base 55.0 g boric acid 9.3 g Na <sub>2</sub> EDTA
TAE (Tris-acetate)	40 mM Tris-acetate (pH 7.6) 1 mM Na <sub>2</sub> EDTA	48.4 Tris-base 11.4 ml glacial acetic acid 20.0 ml 0.5 M Na <sub>2</sub> EDTA (pH 8.0)
TPE (Tris-phosphate)	89 mM Tris-phosphate 2 mM Na <sub>2</sub> EDTA	108.0 g Tris-base 15.5 ml phosphoric acid (85%) 40.0 ml 0.5 M Na <sub>2</sub> EDTA (pH 8.0)

- 0.5 x TBE can be used for agarose gels  
1 x TBE is the standard running buffer for DNA separation in polyacrylamide gels.
- Concentrated stock solutions of TBE tend to develop a precipitate when stored for long periods of time. Store at room temperature and discard any solution that develops a precipitate.

#### Recommended agarose gel concentration for resolving linear DNA

% Agarose	DNA size range (bp)
0.40	2000 – 30,000
0.75	1000 – 15,000
1.00	500 – 10,000
1.25	300 – 5000
1.50	200 – 4000
2.00	100 – 2500

#### Gel loading buffers for non-denaturing agarose gel electrophoresis

(Adapted from Brown, 1991)

Loading buffer	Recipe for 6 x buffer
SBX (store at 15 – 25°C)	40% (w/v) Sucrose 0.25% (w/v) Bromphenol blue 0.25% (w/v) Xylene cyanol
FBX (store at 2 – 8°C)	15% (w/v) Ficoll 400 0.25% (w/v) Bromphenol blue 0.25% (w/v) Xylene cyanol
GBX (store at 2 – 8°C)	30% (v/v) Glycerol 0.25% (w/v) Bromphenol blue 0.25% (w/v) Xylene cyanol

## Denaturing agarose gel electrophoresis

	<b>Alkaline gels</b> (McDonnell et al., 1977)	<b>Glyoxal gels</b> (Bailey and Davidson, 1976)	<b>Formaldehyde gels</b> (Rueger et al., 1996)
<b>Nucleic acid</b>	DNA only	DNA, RNA	RNA (DNA)
<b>Principle</b>	Gels run in alkaline buffer which denatures DNA	Nucleic acids are denatured with glyoxal prior to electrophoresis	Gels are run in the presence of formaldehyde to denature nucleic acids
<b>Gel buffer</b>	50 mM NaCl 1 mM Na <sub>2</sub> EDTA (neutral pH, because at high temperatures alkaline pH hydrolyzes agarose)	10 mM sodium phosphate, pH 7.0 (neutral pH, because glyoxal dissociates from nucleic acids at pH >8.0)	1 x MOPS* 2% (v/v) Formaldehyde
<b>Running buffer</b>	30 mM NaOH, 1 mM Na <sub>2</sub> EDTA (Soak gel in running buffer 30 min before loading.)	Same as gel buffer (The running buffer must be circulated during electrophoresis.)	1 x MOPS*
<b>Loading buffer</b>	Resuspend DNA in: 50 mM NaOH 1 mM Na <sub>2</sub> EDTA 3% (w/v) Ficoll 400 0.025% (w/v) Bromocresol green 0.025% (w/v) Xylene cyanol	Resuspend DNA/RNA in: 8 µl 1 M glyoxal 50% (v/v) DMSO 10 mM sodium phosphate, pH 7.0 Incubate at 50°C for 1 h; cool, then add: 2 µl 10 mM sodium phosphate, pH 7.0 50% (v/v) glycerol 0.4% (w/v) Bromphenol blue	<b>Prepare fresh!!!</b> 250 µl formamide (deionized) 83 µl formaldehyde, 37% (w/v) 50 µl 10 x MOPS-buffer* 0.01% (w/v) Bromphenol blue Add RNA and incubate for 10 min at 65°C; then immediately cool on ice.

\* 10 x MOPS:  
200 mM Morpholinopropanesulfonic acid, pH 7.0  
50 mM sodium acetate  
10 mM Na<sub>2</sub>EDTA

## Nondenaturing polyacrylamide gels: Recommended polyacrylamide gel concentration for resolving dsDNA and behavior of marker dyes

(Adapted from Sambrook et al., 1989)

% Acrylamide	DNA size range (bp)	Bromphenol Blue*	Xylene Cyanol*
3.5	100 – 1000	100	460
5.0	75 – 500	65	260
8.0	50 – 400	45	160
12.0	35 – 250	20	70
15.0	20 – 150	15	60
20.0	5 – 100	12	45

\* Approximate size of DNA fragments (in base pairs) with which the dyes would migrate.

9

## Denaturing polyacrylamide gels: Migration of marker dyes

(Adapted from Sambrook et al., 1989)

% Acrylamide	Bromphenol Blue*	Xylene Cyanol*
5.0	35	130
6.0	26	106
8.0	19	76
10.0	12	55

\* Approximate size of DNA fragments (in nucleotides) with which the dyes would migrate. (Acrylamide: Bisacrylamide (w/w) should be 29:1.)

## Agaroses and molecular weight markers

see inside back cover page.

## Stains and tracker dyes

Compound	MW	Application	Concentration
<b>Acridine Orange</b>	301.8	DNA, RNA stain in gel electrophoresis	1 – 10 µg/ml in staining solution
<b>Bromphenol Blue</b>	670.0	Gel electrophoresis marker dye	0.01 – 0.05% in gel loading buffers
<b>Ethidium bromide</b>	394.3	Fluorescent dye for DNA/ RNA in density gradients and electrophoresis gels	1 µg/ml in running buffer or staining solution
<b>Xylene Cyanol</b>	554.6	Gel electrophoresis marker dye	0.01 – 0.05% in gel loading buffers

## IV. Other Useful Information

### Length of rRNAs from various species

(Lewin, 1987)

Species	RNA	Length (bases)
<i>E. coli</i>	16 S rRNA	1542
	23 S rRNA	2904
Yeast ( <i>S. cerevisiae</i> )	18 S rRNA	2000
	28 S rRNA	3750
<i>Drosophila</i> ( <i>D. melanogaster</i> )	18 S rRNA	2000
	28 S rRNA	4100
Mammal	18 S rRNA	1874
	28 S rRNA	4718