

Human tPA Control Primer Set

(For the amplification of 4.8, 9.3 and 15 kb fragments of the human tPA gene with the Expand Long Template PCR System)

Cat. No. 11 691 104 001

(3 × 25 assays)

Version March 2006

Store at -15 to -25°C

Product description

tPA primers are synthetic oligonucleotide primers designed for the sequence specific amplification of human tPA target DNA (1) by the polymerase chain reaction (PCR).

Concentration 200 ng/μl or 18 pmol/μl in redist. water.

Stability The primers should be stored at -15 to -25°C.

Typical analysis Content from capillary electrophoresis (area %) > 90%. The primers are tested with the Expand¹⁾ Long Template PCR System.

The set contains

- 1. tPA forward 7**
5' primer no. 7, 200 ng/μl (50 μl) for 50 assays
Sequence: 5'-GGA AGT ACA GCT CAG AGT TCT GCA GCA CCC CTG C-3'
- 2. tPA reverse 10**
3' primer no. 10, 200 ng/μl (25 μl) for 25 assays
Sequence: 5'-GAT GCG AAA CTG AGG CTG GCT GTA CTG TCT C-3'
- 3. tPA reverse 14**
3' primer no. 14, 200 ng/μl (25 μl) for 25 assays
Sequence: 5'-CAA AGT CAT GCG GCC ATC GTT CAG ACA CAC C-3'
- 4. tPA forward 1**
5' primer no. 1, 200 ng/μl (25 μl) for 25 assays
Sequence: 5'-CCT TCA CTG TCT GCC TAA CTC CTT CGT GTG TTC C-3'
- 5. tPA reverse 2**
3' primer no. 2, 200 ng/μl (25 μl) for 25 assays
Sequence: 5'-ACT GTG CTT CCT GAC CCA TGG CAG AAG CGC CTT C-3'

Application

The primers are used for the amplification of large PCR fragments from human genomic DNA with the Expand Long Template PCR System*.

Primer concentration for the amplification of tPA fragments is 300 nM (1 μl of the primer in 50 μl reaction).

A detailed protocol is available with the Expand Long Template PCR System* product from Roche Applied Science. Especially for 15 kb amplification the success of PCR depends on the quality of the template used. Therefore, we recommend to use Human Genomic DNA* from Roche Applied Science as control template.

amplified PCR fragments	forward primer gene pos. 5' end	reverse primer gene pos. 5' end	annealing temperature
4.8 kb	tPA forward 7 (26 055)	tPA reverse 10 (30 850)	65°C
9.3 kb	tPA forward 7 (26 055)	tPA reverse 14 (35 388)	65°C
15 kb	tPA forward 1 (8301)	tPA reverse 2 (23 417)	62°C

Reference

- 1 Degen, S. J. et al. (1986) J. Biol. Chem. 261, 6972-6985.

Related products

Product	Pack size	Cat. No.
Expand Long Template PCR System	100 units	11 681 834 001
	500 units	11 681 842 001
Human Genomic DNA	100 μg	11 691 112 001
Expand High Fidelity PCR System	100 units	11 732 641 001
	500 units	11 732 650 001
Human Genomic DNA	100 μg	11 691 112 001

* available from Roche Applied Science

¹⁾ Expand is a trademark of a Member of the Roche Group.

NOTICE TO PURCHASER: DISCLAIMER OF LICENSE

This product is optimized for use in the polymerase chain reaction (PCR) covered by patents owned by F. Hoffmann-La Roche Ltd ("Roche"). No license under these patents to use the PCR Process is conveyed expressly or by implication to the purchaser by the purchase of this product. A license to use the PCR Process for certain research and development activities accompanies the purchase of certain Roche, Applied Biosystems or other licensed suppliers' reagents when used in conjunction with an authorized thermal cycler, or is available from Applied Biosystems. Diagnostic purposes require a license from Roche. Further information on purchasing licenses may be obtained by contacting the Director of Licensing at Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

Contact and Support

To ask questions, solve problems, suggest enhancements or report new applications, please visit our **Online Technical Support Site** at:

www.roche-applied-science.com/support

To call, write, fax, or email us, visit the Roche Applied Science home page, www.roche-applied-science.com, and select your home country. Country-specific contact information will be displayed. Use the Product Search function to find Pack Inserts and Material Safety Data Sheets.



Roche Diagnostics GmbH
Roche Applied Science
68298 Mannheim
Germany