

α -Glucosidase Assay

Colorimetric assay for the quantitative determination of α -glucosidase (maltase) in human semen research samples

Cat. No. 11 742 027 001

30 assays + 10 blanks

Version January 2005

Store at +2 to +8°C

1. What this Product Does

Number of Tests

30 tests and 10 blanks

Kit Contents

Vial	Label	Contents	
1	Reaction Buffer Phosphate buffer, pH 6.8 containing 1% SDS	4 ml	glass bottle blue screw cap
2	Substrate Concentrate Solution containing 4-Nitrophenyl- α -D-glucopyranoside in organic solvent	400 μ l	glass bottle red screw cap
3	Stopping Buffer Sodium carbonate buffer	100 ml	plastic bottle white screw cap
4	Inhibitor Castanospermine, lyophilizate		plastic vessel green screw cap
5	Standard Concentrate Solution containing 4-Nitrophenol in organic solvent	1.5 ml	glass bottle yellow screw cap

Storage and Stability

The kit is stable until expiry date given on the kit label if stored at +2 to +8°C. For stability of reconstituted kit components and working solutions see section 2.1.

Additional Equipment Required

The enzyme incubation during the assay is performed at +37°C in a water bath or other incubator.

To measure the α -glucosidase values a photometer with 1.5 ml cuvettes or a microplate reader with a 405 nm wavelength filter can be used.

A positive displacement pipette (*e.g.*, Gilson Microman) is necessary for pipetting the highly viscous semen samples.

Application

The α -Glucosidase Assay is designed to specifically measure the neutral α -glucosidase in human semen samples **in life science applications. It is not designed for diagnostic applications.** The kit can **not** be used to measure α -glucosidase in other specimens. The assay is compatible with microplate assay format and standard cuvette format.

Assay Time

Approx. 2.5 h

2. How to Use this Product

2.1 Before You Begin

Preparation and Stability of Solutions

Reaction Solution

For each sample to be analyzed dissolve 10 μ l of substrate concentrate solution (bottle 2 warm to +37°C until the content of the bottle is entirely dissolved) in 90 μ l reaction buffer (bottle 1). In addition for each series of samples another 100 μ l reaction solution is necessary for the castanospermine blank. Prewarm the reaction solution to +37°C in a water bath or other incubator before use.

The reaction solution should be always prepared freshly and is stable for several hours.

Inhibitor Solution

Dissolve the lyophilized castanospermine (bottle 4) in 150 μ l redist. water.

The inhibitor solution is stable at +2 to +8°C for two months, if stored separately at -15 to -25°C until the expiry date given on the kit.

Standard Curve

For each series of samples to be analyzed prepare 10 ml of a 100 μ M 4-nitrophenol diluted standard solution by mixing 200 μ l of standard concentrate solution (bottle 5 warm to +37°C until the content of the bottle is entirely dissolved) with 9.8 ml of stopping buffer (bottle 3). Prepare the standard curve according the following table:

μ l 100 μ M 4-nitrophenol (= standard solution)	μ l stopping buffer (bottle 3)	end concentration: 4-nitrophenol standard solution [μ M]
0	1000	0
200	800	20
400	600	40
600	400	60
800	200	80
1000	0	100

2.2 Assay Protocol

- 1 Centrifuge semen sample for 10 min at $1000 \times g$ and transfer the clear supernatant (seminal plasma) to a new tube. The seminal plasma can be used either directly or stored at -15 to -25°C . Repeated freezing and thawing of the seminal plasma does not influence seminal α -glucosidase activity.
- 2 Place $15 \mu\text{l}$ of the semen sample into a 1.5 ml reaction vessel by using a positive displacement pipette (e.g., Gilson Microman). Owing to the high viscosity of the seminal plasma other pipettes will cause higher variances.
Optional: Prepare an internal quality control in the same way (previously measured sample from which aliquots have been stored frozen). Semen samples may be set up in duplicates to decrease variations of the pipetting steps.
- 3 Only to one sample of each series (e.g., to a second sample of the internal quality control pool) add $8 \mu\text{l}$ casanospermine (negative control).
- 4 Add $100 \mu\text{l}$ of prewarmed reaction solution to each sample, mix by vortexing and incubate at $+37^\circ\text{C}$ for 2 h.
- 5 During incubation prepare standard curve.
- 6 Stop incubation by addition of 1 ml stopping buffer and mix by vortexing.
- 7 Transfer entire sample to a 1.5 ml cuvette or a $250 \mu\text{l}$ aliquot into a microplate.
- 8 Read the absorbance of the sample (A_{sample}) and the absorbance of the seminal plasma blank (A_{blank}) respectively at 405 nm against 1 ml stopping buffer. Reading should be completed within 1 h after stopping the reaction.

2.3 Analysis

Calculation

Calculate the net absorbance (ΔA) for every sample:

$$\Delta A = A_{\text{sample}} - A_{\text{blank}}$$

Determine the inverted slope of the standard curve ($\mu\text{M}/\text{absorbance unit}$)

$$V = \text{total volume} = 1.115 \text{ ml}$$

$$v = \text{sample volume} = 0.015 \text{ ml}$$

$$\text{Glucosidase activity in the sample [mU/ml]} = \frac{1}{\text{Slope}} \times \frac{\Delta A \times V}{120 \text{ min} \times v}$$


or:

$$\text{Glucosidase activity in the sample [mU/ml]} = \frac{1}{\text{Slope}} \times \Delta A \times 0,619$$

3. Results

measured sample absorbance (A_{sample})	= 1.45 absorbance units
semen blank absorbance (A_{blank})	= 0.05 absorbance units
Net absorbance of the sample (ΔA)	= 1.4 absorbance units
1/slope of the standard curve	= $60 \mu\text{M}/1.123$ absorbance units = $53.4 \mu\text{M} / \text{absorbance unit}$
Glucosidase activity in the sample	= $53.4 \times 1.4 \times 0.619$ mU/ml = 46.3 mU/ml

⚠ The standard curve should be linear, and so the slope may be determined at any concentration and can be used to multiply all values by the same factor. If the standard curve were to be non-linear and decreasing at higher concentrations (e.g., as may be for some older types of photometers), then a curve should be drawn through the points and the slope is determined directly at the absorbance value of the sample.



4. Additional Information on this Product

4.1 How this Product Works


This product is designed for use in research studies to elucidate the role of α -glucosidase in male infertility. Literature indicates that the measurement of α -glucosidase in human seminal plasma may be used in the evaluation of male infertility. The normal ejaculate is a heterogeneous mixture of the secretions from three organs of the male genital tract. About 50 to 65% of the volume is secreted by the seminal vesicles with about 30 to 40% contributed by the prostate. Secretion of the testes, epididymides and vasa deferentia makes up 3 to 5% of the total ejaculate (1). Fructose (2) and citrate (3) levels in seminal fluid indicate the activity of seminal vesicles and prostate respectively, whereas carnitine (4) and α -glucosidase seem to be correlated with the epididymal function (5). α -Glucosidase activity is easy to measure and seems to give more significant results than carnitine determinations in individuals with proven epididymal dysfunction (6). Although seminal α -glucosidase activity originates mainly from the epididymis the determination of total α -glucosidase activity would also measure α -glucosidase from other organs (e.g., prostate) and thus lead to erroneous results. This test has been substantially improved by taking advantage of the fact that these enzymes have pH-optima in the acidic range and can be totally inhibited by sodium dodecyl sulfate (SDS) (7). A further improvement was achieved by subtracting non- α -glucosidase-regulated degradation of the substrate as a castanospermine inhibited semen blank (8).

Levels of α -glucosidase activity seems to be very low in that cases of azoospermia and asthenozoospermia, where bilateral ductal occlusions are situated between the epididymis and the ejaculatory duct (9, 10). Low levels of α -glucosidase therefore probably reflect a variety of pathological conditions e.g., agenesis of the vas deferens, varicocele, infections, inflammations or vasectomy (9). Normal α -glucosidase activity in subjects with azoospermia and asthenozoospermia indicates according to the literature other causes, such as an arrest of sperm maturation, obstructions between the rete testis and the epididymis or obstructions of the rete testis itself.

The α -Glucosidase Assay is intended as a tool to increase scientific knowledge about the above relationship.

4.2 Test Principle

The α -Glucosidase Assay contains the substrate and all buffers for the measurement of α -glucosidase in human semen or seminal plasma samples. It also contains a standard substance to set up a standard curve. In a first step the semen sample usually is centrifuged and the clear supernatant is used for the assay. Uncleared samples may also be used, but may in some cases lead to higher assay variances. Subsequently the assay is started by adding the substrate in buffer to the sample. The buffer maintains a neutral pH and contains sodium dodecyl sulfate (SDS). Therefore only the neutral α -glucosidase from the epididymis is measured. A blank control sample is set up by mixing the specific inhibitor castanospermine to the sample before adding the substrate. During a 2 hour incubation at +37°C the α -glucosidase containing sample will liberate 4-nitrophenol from the substrate 4-nitrophenyl- α -D-glucopyranoside according to the following equation:



4-Nitrophenol yields a yellow colour upon addition of stopping buffer, which is measured in photometer. After subtraction of the castanospermine inhibited semen blank, the enzyme activity is calculated by comparison to a 4-nitrophenol standard curve.

4.3 Assay Characteristics

Sensitivity

The sensitivity of the assay is 2 mU/ml. Only 15 μ l semen or seminal plasma is necessary due to the high sensitivity.

Measuring Range

The assay is linear between 2 mU/ml and 45 mU/ml.

Precision

Intra-assay variances are only 1–3%, inter-assay variances are up to 5%.

4.4 References

- 1 Wetterauer, U. (1986) *Urol. Res.* **14**, 241.
- 2 Schill, W.-B. (1976) *Med. Klin.* **71**, 1031.
- 3 Hensel, R. & Hornstein, O. P. (1970) *Andrologie* **2**, 61.
- 4 Wetterauer, U. & Heite, H.-J. (1976) *Akt. Dermatol.* **3**, 93.
- 5 Paquin, R. et al. (1984) *J. Androl.* **5**, 277.
- 6 Fourie, M. H. et al. (1993) *Arch. Androl.* **31**, 209.
- 7 Tremblay, R. R. (1986) In: *Andrology, Male Fertility and Sterility* (Paulson, J. D. et al. eds.) Academic Press, London, p. 273.
- 8 Cooper, T. G. et al. (1990) *Int. J. Androl.* **13**, 297.
- 9 Cooper, T. G. et al. (1988) *J. Androl.* **9**, 91.
- 10 Viljoen, M. H. et al. (1990) *Andrologia* **22**, 205.

5. Supplementary Information

5.1 Conventions


Text Conventions

To make information consistent and memorable, the following text conventions are used in this package insert:

Text Convention	Use
Numbered Instructions labeled ① ②, etc.	Steps in a procedure that must be performed in the order listed
Asterisk *	Denotes a product available from Roche Applied Science

Symbols

In this package insert the following symbol is used to highlight important information:

Symbol	Description
	Important Note: Information critical to the success of the procedure or use of the product.

5.2 Changes to Previous Version

Version January 2005 is completely revised.

Contact and Support

To ask questions, solve problems, suggest enhancements or report new applications, please visit our **Online Technical Support** Site at:

www.roche-applied-science.com/support

To call, write, fax, or email us, visit the Roche Applied Science home page, www.roche-applied-science.com, and select your home country.

Country-specific contact information will be displayed.

Use the Product Search function to find Pack Inserts and Material Safety Data Sheets.

Inspiring Discovery



Diagnostics

Roche Diagnostics GmbH
Roche Applied Science
68298 Mannheim
Germany