

## The DIG System

*Labeling and Detection of Nucleic Acids*




## The DIG System

### Specifically Label and Detect Nucleic Acids

*Publishable results require high level specific detection and low background. Do your hybridizations have nonspecific signals and high background? The DIG System is ideal for nucleic acid labeling. Flexibly use colorimetric, luminescent or fluorescent signal detection. Achieve high sensitivity and low background in very short exposure times*

- **Specificity:** DIG antibodies do not bind other substrates.
- **Versatility:** Use DIG labeled probes for filters and *in situ* hybridization.
- **Proven:** Thousands of publications show why DIG is superior to radioactivity.



**Figure 1: DIG-dUTP coupling is via an alkali-labile ester bond.** Anti-DIG antibodies only bind to DIG, not other biological substrates. The alkali-labile ester bond ideal for filter stripping and reprobing.

Products are for life science research only. Not for use in diagnostic procedures.

#### Trademarks

DIG EASY HYB and TELOTAGGG are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.


#### License Disclaimer

Manufactured under license from Geron Corporation.

## Combine DIG Labeling with Established Protocols

*High specificity and sensitivity are the reason, researchers worldwide choose the DIG System to detect nucleic acids using filter and in situ hybridization. Rely on the DIG System for straightforward, sensitive, and efficient nonradioactive labeling and detection. Use robust procedures and established protocols for low background and a high signal-to-noise ratio.*

Probes for NA-labeling	Methods for Labeling	Hybridization and Detection
DIG labeling produces high sensitivity and low background.  Easily label <ul style="list-style-type: none"> <li>- DNA</li> <li>- RNA</li> <li>- oligonucleotides</li> </ul>	Label the probe using familiar, efficient procedures.  Quickly achieve superb results using PCR labeling or <i>in vitro</i> transcription.	The DIG System can be applied to filter or <i>in situ</i> hybridization.  Established easy-to-use protocols make stripping and reprobing straightforward.



**Figure 2: Example detecting DIG-labeled nucleic acids using chemiluminescence substrates.**

Immobilized target nucleic acids are hybridized to a DIG-labeled probe. Subsequent detection is performed using high affinity anti-digoxigenin antibodies coupled either to alkaline phosphatase (AP), horseradish peroxidase (POD), and fluorescein or rhodamine for colorimetric, chemiluminescence or fluorescence detection.

# The DIG System

## Find Applications and Corresponding Products


### ■ Blotting Applications

Simple, effective nonradioactive labeling and detection of nucleic acids ensures high signal-to-noise. All Kits are QC (quality control) tested, to function in different blot applications and guaranteed to be DNase and RNase free

according to the current quality procedures. With DIG-labeled probes, you can easily detect single-copy genes on Southern blots, unique mRNAs on northern blots, or rare recombinants in bacterial colonies or viral plaques.

Labeling	Immobilization	Hybridization and Detection
<b>PCR</b> <ul style="list-style-type: none"> <li>■ DIG Probe Synthesis Kit</li> <li>■ DIG Northern Starter Kit</li> </ul> <b>In Vitro Transcription</b> <ul style="list-style-type: none"> <li>■ DIG RNA Labeling Kit (SP6/T7)</li> <li>■ DIG RNA Labeling Mix</li> </ul>	<ul style="list-style-type: none"> <li>■ Nylon Membranes, positively charged</li> <li>■ Buffers in a Box</li> <li>■ Molecular Weight Marker, DIG-labeled (DNA or RNA)</li> </ul>	<ul style="list-style-type: none"> <li>■ DIG Easy Hyb</li> <li>■ Hybridization Bags</li> <li>■ Actin RNA Probe Labeled DIG (as control)</li> <li>■ Anti-Digoxigenin-AP, Fab fragments</li> <li>■ CDP-<i>Star</i>, ready-to-use</li> <li>■ CDP-<i>Star</i>, ready-to-use NBT/BCIP</li> <li>■ DIG Wash and Block Buffer Set</li> <li>■ Lumi Film</li> </ul>

■ Products are recommended but not absolutely required.



**Figure 3: Northern blot hybridized with a DIG-labeled p21 probe.**

A strong increase over time is visible for the p21 mRNA. Exposure time: 3 minutes after 2 hours incubation with CDP-*Star*.

#### → Find additional information

Obtain Cancer Research Application Note No. 10,  
DIG Application Manual for Filter Hybridization,  
Catalog Number 05 353 149 001


## ■ *In Situ* Applications

*In situ* hybridization techniques allow specific nucleic acid sequences to be detected in morphologically preserved chromosomes, cells, or tissue sections. In combination with immunocytochemistry, *in situ* hybridization can relate microscopic topological

information to gene activity at the DNA, mRNA, and protein levels. The DIG System can be used for fresh-frozen or fixed tissue and cells to detect expression patterns or chromosome aberrations.

Labeling	Immobilization	Hybridization and Detection
<p><b>Nick Translation</b></p> <ul style="list-style-type: none"> <li>■ Nick Translation Mix</li> <li>■ DIG-Nick Translation Mix</li> </ul> <p><b>PCR</b></p> <ul style="list-style-type: none"> <li>■ PCR DIG Probe Synthesis Kit</li> </ul> <p><b><i>In Vitro</i> Transcription</b></p> <ul style="list-style-type: none"> <li>■ DIG RNA Labeling Kit (SP6/T7)</li> <li>■ DIG RNA Labeling Mix</li> </ul>	<p><b>Tissue Preparation</b> depending on your laboratory application</p> <ul style="list-style-type: none"> <li>- FFPE</li> <li>- Fresh Frozen</li> <li>- Whole Mounts</li> <li>- Chromosome Spreads</li> </ul>	<ul style="list-style-type: none"> <li>■ DIG Wash and Block Buffer Set</li> <li>■ Anti DIG AP Antibody, Fab fragments</li> <li>■ Anti-Digoxigenin-Fluorescein, Fab fragments</li> <li>■ Anti-Digoxigenin-Rhodamine, Fab fragments</li> <li>■ Anti-Digoxigenin-POD, Fab fragments</li> <li>■ NBT/BCIP/HNPP</li> <li>■ POD conjugates</li> </ul>

The DIG System can also be used in sophisticated micro-RNA staining techniques using common model organisms.



**Figure 4: Dorsal view of a zebrafish embryo whole mount *in situ* hybridization.** The slide was stained with a DIG-labeled RNA probe for the proneural transcription factor *atoh1b* and detected with NBT/BCIP.

→ [Find additional information](#) in the DIG Application Manual for Nonradioactive *In Situ* Hybridization, 4<sup>th</sup> Edition, Catalog Number 05 353 122 001

## ■ Additional Applications

Following nonradioactive techniques go beyond standard nucleic acid labeling and detection. Highly efficient and sensitive techniques include:

- Glycan labeling
- Protein labeling
- Gel Shift assay
- PCR ELISA
- Telomere assay

→ [Find additional information](#) about nonradioactive techniques for protein labeling, western blotting, and other applications, [www.roche-applied-science.com](http://www.roche-applied-science.com)

## Ordering Information

	<b>Catalog Number</b>	<b>Pack Size</b>
<b>Kits for Labeling and Detection</b>		
<b>DIG-High Prime DNA Labeling and Detection Starter Kit I</b>	11 745 832 910	1 kit (12 labeling reactions and 24 detection reactions)
<b>DIG-High Prime DNA Labeling and Detection Starter Kit II</b>	11 585 614 910	1 kit (12 labeling reactions and 24 detection reactions)
<b>DIG DNA Labeling and Detection Kit</b>	11 093 657 910	1 kit (25 labeling reactions and 50 blots)
<b>DIG Northern Starter Kit</b>	12 039 672 910	1 kit (10 labeling reactions and detection of 10 blots 10 x 10 cm <sup>2</sup> )
<b>Kits for Labeling</b>		
<b>PCR DIG Probe Synthesis Kit</b>	11 636 090 910	1 kit (25 labeling reactions)
<b>DIG DNA Labeling Kit</b>	11 175 033 910	1 kit (40 labeling reactions)
<b>DIG-High Prime</b>	11 585 606 910	160 µl (40 labeling reactions)
<b>DIG RNA Labeling Kit (SP6/T7)</b>	11 175 025 910	1 kit (2 x 10 labeling reactions)
<b>DIG Oligonucleotide 3'-End Labeling Kit, 2<sup>nd</sup> generation</b>	03 353 575 910	1 kit (25 labeling reactions)
<b>DIG Oligonucleotide 5'-End Labeling Set</b>	11 480 863 001	1 set (10 labeling reactions)
<b>DIG Oligonucleotide Tailing Kit, 2<sup>nd</sup> generation</b>	03 353 583 910	1 kit (25 tailing reactions)
<b>Mixes for Labeling</b>		
<b>DIG DNA Labeling Mix</b>	11 277 065 910	50 µl (25 reactions)
<b>DIG-High Prime</b>	11 585 606 910	160 µl (40 reactions)
<b>PCR DIG Labeling Mix**</b>	11 585 550 910	2 x 250 µl (2 x 25 reactions)
<b>DIG RNA Labeling Mix</b>	11 277 073 910	40 µl (20 reactions)
<b>Biotin-High Prime</b>	11 585 649 910	100 µl (25 reactions)
<b>Biotin RNA Labeling Mix</b>	11 685 597 910	40 µl (20 reactions)
<b>Fluorescein-High Prime*</b>	11 585 622 910	100 µl (25 reactions)
<b>PCR Fluorescein Labeling Mix</b>	11 636 154 910	10 PCR reactions (100 µl each)
<b>Fluorescein RNA Labeling Mix</b>	11 685 619 910	40 µl (20 reactions)
<b>Nick Translation Mix*</b>	11 745 808 910	200 µl
<b>DIG-Nick Translation Mix*</b>	11 745 816 910	160 µl
<b>Biotin-Nick Translation Mix*</b>	11 745 824 910	160 µl
<b>Nucleotides for Labeling</b>		
<b>Digoxigenin-11-dUTP, alkali-labile</b>	11 573 152 910	25 nmol (25 µl)
	11 573 179 910	125 nmol (125 µl)
	11 093 088 910	25 nmol (25 µl)
	11 558 706 910	125 nmol (125 µl)
<b>Digoxigenin-11-dUTP, alkali-stable</b>	11 570 013 910	5 x 125 nmol (5 x 125 µl)
	11 209 256 910	250 nmol (25 µl)
	03 359 247 910	200 nmol (57 µl)
<b>Digoxigenin-11-UTP</b>	11 363 905 910	25 nmol (25 µl)
<b>Digoxigenin-11-ddUTP</b>	11 363 905 910	25 nmol (25 µl)
<b>Biotin-16-dUTP</b>	11 093 070 910	50 nmol (50 µl)
<b>Biotin-16-UTP</b>	11 388 908 910	250 nmol (25 µl)
<b>Biotin-16-ddUTP</b>	11 427 598 910	25 nmol (25 µl)
<b>Fluorescein-12-dUTP*</b>	11 373 242 910	25 nmol (25 µl)
<b>Fluorescein-12-UTP*</b>	11 427 857 910	250 nmol (25 µl)
<b>Tetramethyl-rhodamine-5-dUTP*</b>	11 534 378 910	25 nmol (25 µl)
<b>Kits for Detection</b>		
<b>DIG Luminescent Detection Kit</b>	11 363 514 910	1 kit (50 blots)
<b>DIG Nucleic Acid Detection Kit</b>	11 175 041 910	1 kit (40 blots)
<b>Anti-DIG Antibody for Detection on Membranes</b>		
<b>Anti-Digoxigenin-AP, Fab fragments</b>	11 093 274 910	150 U (200 µl)
<b>Anti-Digoxigenin-Fluorescein, Fab fragments*</b>	11 207 741 910	200 µg
<b>Anti-Digoxigenin-Rhodamine, Fab fragments*</b>	11 207 750 910	200 µg
<b>Anti-Digoxigenin-POD, Fab fragments</b>	11 207 733 910	150 U

Chemiluminescent Substrates for Alkaline Phosphatase	Catalog Number	Pack Size
<b>CDP-Star</b>	11 685 627 001	1 ml
	11 759 051 001	2 x 1 ml
<b>CDP-Star, ready-to-use</b>	12 041 677 001	2 x 50 ml
<b>CSPD</b>	11 655 884 001	1 ml
<b>CSPD, ready-to-use</b>	11 755 633 001	2 x 50 ml
Chromogenic Substrates for Alkaline Phosphatase		
<b>BM purple*</b>	11 442 074 001	100 ml
<b>5-Bromo-4-chloro-3-indolyl-phosphate (BCIP)*</b>	11 383 221 001	3 ml (150 mg)
<b>HNPP Fluorescent Detection Set*</b>	11 758 888 001	1 set (5 mg HNPP, 100 mg Fast Red TR)
<b>NBT (4-Nitroblue tetrazolium chloride)*</b>	11 383 213 001	3 ml (300 mg)
<b>NBT/BCIP Ready-to-Use Tablets*</b>	11 697 471 001	20 tablets
<b>NBT/BCIP Stock Solution*</b>	11 681 451 001	8 ml
<b>Fast Red Tablets*</b>	11 496 549 001	20 tablets
Molecular Weight Markers, DIG-labeled		
<b>DNA Molecular Weight Marker II, DIG-labeled</b>	11 218 590 910	5 µg (500 µl)
<b>DNA Molecular Weight Marker III, DIG-labeled</b>	11 218 603 910	5 µg (500 µl)
<b>DNA Molecular Weight Marker VI, DIG-labeled</b>	11 218 611 910	5 µg (500 µl)
<b>DNA Molecular Weight Marker VII, DIG-labeled</b>	11 669 940 910	5 µg (500 µl)
<b>DNA Molecular Weight Marker VIII, DIG-labeled</b>	11 449 451 910	5 µg (500 µl)
<b>RNA Molecular Weight Marker I, DIG-labeled</b>	11 526 529 910	4 µg (200 µl)
<b>RNA Molecular Weight Marker II, DIG-labeled</b>	11 526 537 910	2 µg (200 µl)
<b>RNA Molecular Weight Marker III, DIG-labeled</b>	11 373 099 910	2 µg (200 µl)
Nucleic Acids and Probes, DIG-labeled		
<b>Actin RNA Probe, DIG-labeled</b>	11 498 045 910	2 µg
<b>DIG-labeled Control DNA</b>	11 585 738 910	50 µl
<b>DIG-labeled Control RNA</b>	11 585 746 910	50 µl
Additional Kits for Nonradioactive Analysis		
<b>PCR ELISA, DIG Detection</b>	11 636 111 910	192 detection reactions
<b>DIG Gel Shift Kit, 2nd generation</b>	03 353 591 910	1 kit
<b>Lumi-Light<sup>PLUS</sup> Western Blotting Kit</b>	12 015 218 001	1 kit
<b>TeloTAGGG Telomere Length Assay</b>	12 209 136 001	1 kit
Membranes and Films for Nonradioactive Analysis		
<b>Nylon Membranes, positively charged</b>	11 209 272 001	10 sheets (20 x 30 cm)
	11 209 299 001	20 sheets (10 x 15 cm)
	11 417 240 001	1 roll (0.3 x 3 m)
<b>Nylon Membranes for Colony and Plaque Hybridization</b>	11 699 075 001	50 discs (each 82 mm diameter)
	11 699 083 001	50 discs (each 132 mm diameter)
<b>Lumi-Film Chemiluminescent Detection Film</b>	11 666 916 001	100 films (7.1 x 9.4 inches 18 x 24 cm)
	11 666 657 001	100 films (8 x 10 inches 20.3 x 25.4 cm)
<b>Hybridization Bags</b>	11 666 649 001	50 bags
Additional Reagents for Nonradioactive Analysis		
<b>DIG Wash and Block Buffer Set</b>	11 585 762 001	1 set (30 blots)
<b>DIG Easy Hyb</b>	11 603 558 001	500 ml
<b>DIG Easy Hyb Granules</b>	11 796 895 001	Granules for 6 x 100 ml
<b>Buffers in a Box, SSC, 20x</b>	11 666 681 001	4 l
<b>Blocking Reagent</b>	11 096 176 001	50 g
<b>DNA, MB-grade (fish sperm)*</b>	11 467 140 001	500 mg (50 ml)
<b>COT Human DNA</b>	11 581 074 001	500 µg (500 µl)
<b>tRNA, RNase negative*</b>	10 109 541 001	100 mg
	10 109 550 001	500 mg
<b>RNA, from yeast*</b>	10 109 223 001	100 g
<b>Hexanucleotide Mix</b>	11 277 081 001	100 µl (50 reactions)
<b>Digoxigenin-3-O-methylcarbonyl-ε-aminocaproic-acid-N-hydroxy-succinimide ester</b>	11 333 054 001	5 mg
<b>DAB Substrate*</b>	11 718 096 001	1 pack

\* Recommended for *in situ* hybridization.

\*\* Recommended for direct detection.

For ordering information visit our homepage: [www.roche-applied-science.com](http://www.roche-applied-science.com)

For life science research only. Not for use in diagnostic procedures.

**Published by**

Roche Diagnostics GmbH  
Sandhofer Straße 116  
68305 Mannheim  
Germany

© 2011 Roche Diagnostics.  
All rights reserved.

[www.roche-applied-science.com](http://www.roche-applied-science.com)

① 1211

