

cOmplete Protease Inhibitor Cocktail Tablets

Now Available in Even More Convenient *EASYPacks*

Protein protection against protease activity is an important step during isolation or purification of proteins. For complete protection against a multitude of proteases different protease inhibitors are needed. The time-consuming and cumbersome job of weighing small amounts of different protease inhibitors on an analytical scale, and dissolving the mix in DMSO is avoided by using cOmplete Protease Inhibitor Cocktail Tablets.



cOmplete *EASYPack* foil strips make cOmplete Tablets even more convenient than ever before

cOmplete is a proprietary blend of protease inhibitors, formulated as a ready-to-use water soluble tablet. After addition of the tablet to the lysis buffer, the proteins can be instantly protected against a broad range of proteases, which eliminates the time-consuming search for the right protease inhibitors.

Protein protection using cOmplete is now even more convenient with the new *EASYPacks*. Each tablet is individually packed in a foil strip and can be used for:

- ➔ Consistent inhibition of a multitude of protease classes, including serine proteases, cysteine proteases, and metalloproteases
- ➔ Inhibition of proteolytic activity in extracts from almost any tissue or cell type, including animals, plants, yeast, bacteria, and fungi.
For tissue and cell type examples visit www.roche-applied-science.com/proteaseinhibitor
- ➔ Delivery of a consistent dose of protease inhibitors every time
- ➔ Stable, nontoxic protection in aqueous buffers

- ➔ Maintenance of the stability of metal-dependent proteins, and function of purification techniques (*i.e.*, IMAC [immobilized metal-affinity chromatography] for isolation of Poly-His-tagged proteins) by using EDTA-free cOmplete Protease Inhibitor Cocktail Tablets

cOmplete tablets are available in glass vials and *EASYPacks* in two tablet sizes for either 10 ml or 50 ml of lysate, with or without EDTA.

Ask your local Roche sales representative for a two-tablets trial pack of cOmplete, Mini or cOmplete, Mini, EDTA-free and see how convenient and reliable protease inhibition can be.



For more information on cOmplete and our other products for protease inhibition, visit

www.roche-applied-science.com/proteaseinhibitor



cOmplete Protease Inhibitor Cocktail Tablets supplied in *EASYPacks*:

Product	Pack Size	Cat. No.
cOmplete	20 tablets (each for 50 ml)	04 693 116 001
cOmplete, Mini	30 tablets (each for 10 ml)	04 693 124 001
cOmplete, EDTA-free	20 tablets (each for 50 ml)	04 693 132 001
cOmplete, Mini, EDTA-free	30 tablets (each for 10 ml)	04 693 159 001